

SHOCKWAVE™

Kategorivelger i forhold til brukervekt og aktivitetsnivå

Aktivitet / Vekt	44 - 52 kg	53 - 59 kg	60 - 68 kg	69 - 77 kg	78 - 88 kg	89 - 100 kg	101 - 116 kg	117 - 130 kg	131 - 147 kg
Lavt	1	1	2	3	4	5	6	7	8
Normalt	1	2	3	4	5	6	7	8	-
Høyt	2	3	4	5	6	7	8	-	-

Art.nr.	Modell	Lev.art.nr.
722402310	SHOCKWAVE FOOT, CAT.1 SIZE 23CM RU	F20-S3-01A23-RU
722402311	SHOCKWAVE FOOT, CAT.1 SIZE 23CM SL	F20-S3-01A23-SL
722402312	SHOCKWAVE FOOT, CAT.1 SIZE 23CM SR	F20-S3-01A23-SR
722402410	SHOCKWAVE FOOT, CAT.1 SIZE 24CM RU	F20-S3-01A24-RU
722402411	SHOCKWAVE FOOT, CAT.1 SIZE 24CM SL	F20-S3-01A24-SL
722402412	SHOCKWAVE FOOT, CAT.1 SIZE 24CM SR	F20-S3-01A24-SR
722402510	SHOCKWAVE FOOT, CAT.1 SIZE 25CM RU	F20-S3-01A25-RU
722402511	SHOCKWAVE FOOT, CAT.1 SIZE 25CM SL	F20-S3-01A25-SL
722402512	SHOCKWAVE FOOT, CAT.1 SIZE 25CM SR	F20-S3-01A25-SR
722402610	SHOCKWAVE FOOT, CAT.1 SIZE 26CM RU	F20-S3-01A26-RU
722402611	SHOCKWAVE FOOT, CAT.1 SIZE 26CM SL	F20-S3-01A26-SL
722402612	SHOCKWAVE FOOT, CAT.1 SIZE 26CM SR	F20-S3-01A26-SR
722402710	SHOCKWAVE FOOT, CAT.1 SIZE 27CM RU	F20-S3-01A27-RU
722402711	SHOCKWAVE FOOT, CAT.1 SIZE 27CM SL	F20-S3-01A27-SL
722402712	SHOCKWAVE FOOT, CAT.1 SIZE 27CM SR	F20-S3-01A27-SR
722402810	SHOCKWAVE FOOT, CAT.1 SIZE 28CM RU	F20-S3-01A28-RU
722402811	SHOCKWAVE FOOT, CAT.1 SIZE 28CM SL	F20-S3-01A28-SL
722402812	SHOCKWAVE FOOT, CAT.1 SIZE 28CM SR	F20-S3-01A28-SR
722402910	SHOCKWAVE FOOT, CAT.1 SIZE 29CM RU	F20-S3-01A29-RU
722403010	SHOCKWAVE FOOT, CAT.1 SIZE 30CM RU	F20-S3-01A30-RU
722402320	SHOCKWAVE FOOT, CAT.2 SIZE 23CM RU	F20-S3-02A23-RU
722402321	SHOCKWAVE FOOT, CAT.2 SIZE 23CM SL	F20-S3-02A23-SL
722402322	SHOCKWAVE FOOT, CAT.2 SIZE 23CM SR	F20-S3-02A23-SR
722402420	SHOCKWAVE FOOT, CAT.2 SIZE 24CM RU	F20-S3-02A24-RU
722402421	SHOCKWAVE FOOT, CAT.2 SIZE 24CM SL	F20-S3-02A24-SL
722402422	SHOCKWAVE FOOT, CAT.2 SIZE 24CM SR	F20-S3-02A24-SR
722402520	SHOCKWAVE FOOT, CAT.2 SIZE 25CM RU	F20-S3-02A25-RU
722402521	SHOCKWAVE FOOT, CAT.2 SIZE 25CM SL	F20-S3-02A25-SL
722402522	SHOCKWAVE FOOT, CAT.2 SIZE 25CM SR	F20-S3-02A25-SR
722402620	SHOCKWAVE FOOT, CAT.2 SIZE 26CM RU	F20-S3-02A26-RU
722402621	SHOCKWAVE FOOT, CAT.2 SIZE 26CM SL	F20-S3-02A26-SL
722402622	SHOCKWAVE FOOT, CAT.2 SIZE 26CM SR	F20-S3-02A26-SR
722402720	SHOCKWAVE FOOT, CAT.2 SIZE 27CM RU	F20-S3-02A27-RU
722402721	SHOCKWAVE FOOT, CAT.2 SIZE 27CM SL	F20-S3-02A27-SL
722402722	SHOCKWAVE FOOT, CAT.2 SIZE 27CM SR	F20-S3-02A27-SR
722402820	SHOCKWAVE FOOT, CAT.2 SIZE 28CM RU	F20-S3-02A28-RU
722402821	SHOCKWAVE FOOT, CAT.2 SIZE 28CM SL	F20-S3-02A28-SL
722402822	SHOCKWAVE FOOT, CAT.2 SIZE 28CM SR	F20-S3-02A28-SR

722402920	SHOCKWAVE FOOT, CAT.2 SIZE 29CM RU	F20-S3-02A29-RU
722403020	SHOCKWAVE FOOT, CAT.2 SIZE 30CM RU	F20-S3-02A30-RU
722402330	SHOCKWAVE FOOT, CAT.3 SIZE 23CM RU	F20-S3-03A23-RU
722402331	SHOCKWAVE FOOT, CAT.3 SIZE 23CM SL	F20-S3-03A23-SL
722402332	SHOCKWAVE FOOT, CAT.3 SIZE 23CM SR	F20-S3-03A23-SR
722402430	SHOCKWAVE FOOT, CAT.3 SIZE 24CM RU	F20-S3-03A24-RU
722402431	SHOCKWAVE FOOT, CAT.3 SIZE 24CM SL	F20-S3-03A24-SL
722402432	SHOCKWAVE FOOT, CAT.3 SIZE 24CM SR	F20-S3-03A24-SR
722402530	SHOCKWAVE FOOT, CAT.3 SIZE 25CM RU	F20-S3-03A25-RU
722402531	SHOCKWAVE FOOT, CAT.3 SIZE 25CM SL	F20-S3-03A25-SL
722402532	SHOCKWAVE FOOT, CAT.3 SIZE 25CM SR	F20-S3-03A25-SR
722402630	SHOCKWAVE FOOT, CAT.3 SIZE 26CM RU	F20-S3-03A26-RU
722402631	SHOCKWAVE FOOT, CAT.3 SIZE 26CM SL	F20-S3-03A26-SL
722402632	SHOCKWAVE FOOT, CAT.3 SIZE 26CM SR	F20-S3-03A26-SR
722402730	SHOCKWAVE FOOT, CAT.3 SIZE 27CM RU	F20-S3-03A27-RU
722402731	SHOCKWAVE FOOT, CAT.3 SIZE 27CM SL	F20-S3-03A27-SL
722402732	SHOCKWAVE FOOT, CAT.3 SIZE 27CM SR	F20-S3-03A27-SR
722402830	SHOCKWAVE FOOT, CAT.3 SIZE 28CM RU	F20-S3-03A28-RU
722402831	SHOCKWAVE FOOT, CAT.3 SIZE 28CM SL	F20-S3-03A28-SL
722402832	SHOCKWAVE FOOT, CAT.3 SIZE 28CM SR	F20-S3-03A28-SR
722402930	SHOCKWAVE FOOT, CAT.3 SIZE 29CM RU	F20-S3-03A29-RU
722403030	SHOCKWAVE FOOT, CAT.3 SIZE 30CM RU	F20-S3-03A30-RU
722402340	SHOCKWAVE FOOT, CAT.4 SIZE 23CM RU	F20-S3-04A23-RU
722402341	SHOCKWAVE FOOT, CAT.4 SIZE 23CM SL	F20-S3-04A23-SL
722402342	SHOCKWAVE FOOT, CAT.4 SIZE 23CM SR	F20-S3-04A23-SR
722402440	SHOCKWAVE FOOT, CAT.4 SIZE 24CM RU	F20-S3-04A24-RU
722402441	SHOCKWAVE FOOT, CAT.4 SIZE 24CM SL	F20-S3-04A24-SL
722402442	SHOCKWAVE FOOT, CAT.4 SIZE 24CM SR	F20-S3-04A24-SR
722402540	SHOCKWAVE FOOT, CAT.4 SIZE 25CM RU	F20-S3-04A25-RU
722402541	SHOCKWAVE FOOT, CAT.4 SIZE 25CM SL	F20-S3-04A25-SL
722402542	SHOCKWAVE FOOT, CAT.4 SIZE 25CM SR	F20-S3-04A25-SR
722402640	SHOCKWAVE FOOT, CAT.4 SIZE 26CM RU	F20-S3-04A26-RU
722402641	SHOCKWAVE FOOT, CAT.4 SIZE 26CM SL	F20-S3-04A26-SL
722402642	SHOCKWAVE FOOT, CAT.4 SIZE 26CM SR	F20-S3-04A26-SR
722402740	SHOCKWAVE FOOT, CAT.4 SIZE 27CM RU	F20-S3-04A27-RU
722402741	SHOCKWAVE FOOT, CAT.4 SIZE 27CM SL	F20-S3-04A27-SL
722402742	SHOCKWAVE FOOT, CAT.4 SIZE 27CM SR	F20-S3-04A27-SR
722402840	SHOCKWAVE FOOT, CAT.4 SIZE 28CM RU	F20-S3-04A28-RU
722402841	SHOCKWAVE FOOT, CAT.4 SIZE 28CM SL	F20-S3-04A28-SL
722402842	SHOCKWAVE FOOT, CAT.4 SIZE 28CM SR	F20-S3-04A28-SR
722402940	SHOCKWAVE FOOT, CAT.4 SIZE 29CM RU	F20-S3-04A29-RU
722403040	SHOCKWAVE FOOT, CAT.4 SIZE 30CM RU	F20-S3-04A30-RU
722402350	SHOCKWAVE FOOT, CAT.5 SIZE 23CM RU	F20-S3-05A23-RU
722402351	SHOCKWAVE FOOT, CAT.5 SIZE 23CM SL	F20-S3-05A23-SL
722402352	SHOCKWAVE FOOT, CAT.5 SIZE 23CM SR	F20-S3-05A23-SR
722402450	SHOCKWAVE FOOT, CAT.5 SIZE 24CM RU	F20-S3-05A24-RU
722402451	SHOCKWAVE FOOT, CAT.5 SIZE 24CM SL	F20-S3-05A24-SL
722402452	SHOCKWAVE FOOT, CAT.5 SIZE 24CM SR	F20-S3-05A24-SR

722402550	SHOCKWAVE FOOT, CAT.5 SIZE 25CM RU	F20-S3-05A25-RU
722402551	SHOCKWAVE FOOT, CAT.5 SIZE 25CM SL	F20-S3-05A25-SL
722402552	SHOCKWAVE FOOT, CAT.5 SIZE 25CM SR	F20-S3-05A25-SR
722402650	SHOCKWAVE FOOT, CAT.5 SIZE 26CM RU	F20-S3-05A26-RU
722402651	SHOCKWAVE FOOT, CAT.5 SIZE 26CM SL	F20-S3-05A26-SL
722402652	SHOCKWAVE FOOT, CAT.5 SIZE 26CM SR	F20-S3-05A26-SR
722402750	SHOCKWAVE FOOT, CAT.5 SIZE 27CM RU	F20-S3-05A27-RU
722402751	SHOCKWAVE FOOT, CAT.5 SIZE 27CM SL	F20-S3-05A27-SL
722402752	SHOCKWAVE FOOT, CAT.5 SIZE 27CM SR	F20-S3-05A27-SR
722402850	SHOCKWAVE FOOT, CAT.5 SIZE 28CM RU	F20-S3-05A28-RU
722402851	SHOCKWAVE FOOT, CAT.5 SIZE 28CM SL	F20-S3-05A28-SL
722402852	SHOCKWAVE FOOT, CAT.5 SIZE 28CM SR	F20-S3-05A28-SR
722402950	SHOCKWAVE FOOT, CAT.5 SIZE 29CM RU	F20-S3-05A29-RU
722403050	SHOCKWAVE FOOT, CAT.5 SIZE 30CM RU	F20-S3-05A30-RU
722402360	SHOCKWAVE FOOT, CAT.6 SIZE 23CM RU	F20-S3-06A23-RU
722402361	SHOCKWAVE FOOT, CAT.6 SIZE 23CM SL	F20-S3-06A23-SL
722402362	SHOCKWAVE FOOT, CAT.6 SIZE 23CM SR	F20-S3-06A23-SR
722402460	SHOCKWAVE FOOT, CAT.6 SIZE 24CM RU	F20-S3-06A24-RU
722402461	SHOCKWAVE FOOT, CAT.6 SIZE 24CM SL	F20-S3-06A24-SL
722402462	SHOCKWAVE FOOT, CAT.6 SIZE 24CM SR	F20-S3-06A24-SR
722402560	SHOCKWAVE FOOT, CAT.6 SIZE 25CM RU	F20-S3-06A25-RU
722402561	SHOCKWAVE FOOT, CAT.6 SIZE 25CM SL	F20-S3-06A25-SL
722402562	SHOCKWAVE FOOT, CAT.6 SIZE 25CM SR	F20-S3-06A25-SR
722402660	SHOCKWAVE FOOT, CAT.6 SIZE 26CM RU	F20-S3-06A26-RU
722402661	SHOCKWAVE FOOT, CAT.6 SIZE 26CM SL	F20-S3-06A26-SL
722402662	SHOCKWAVE FOOT, CAT.6 SIZE 26CM SR	F20-S3-06A26-SR
722402760	SHOCKWAVE FOOT, CAT.6 SIZE 27CM RU	F20-S3-06A27-RU
722402761	SHOCKWAVE FOOT, CAT.6 SIZE 27CM SL	F20-S3-06A27-SL
722402762	SHOCKWAVE FOOT, CAT.6 SIZE 27CM SR	F20-S3-06A27-SR
722402860	SHOCKWAVE FOOT, CAT.6 SIZE 28CM RU	F20-S3-06A28-RU
722402861	SHOCKWAVE FOOT, CAT.6 SIZE 28CM SL	F20-S3-06A28-SL
722402862	SHOCKWAVE FOOT, CAT.6 SIZE 28CM SR	F20-S3-06A28-SR
722402960	SHOCKWAVE FOOT, CAT.6 SIZE 29CM RU	F20-S3-06A29-RU
722403060	SHOCKWAVE FOOT, CAT.6 SIZE 30CM RU	F20-S3-06A30-RU
722402370	SHOCKWAVE FOOT, CAT.7 SIZE 23CM RU	F20-S3-07A23-RU
722402371	SHOCKWAVE FOOT, CAT.7 SIZE 23CM SL	F20-S3-07A23-SL
722402372	SHOCKWAVE FOOT, CAT.7 SIZE 23CM SR	F20-S3-07A23-SR
722402470	SHOCKWAVE FOOT, CAT.7 SIZE 24CM RU	F20-S3-07A24-RU
722402471	SHOCKWAVE FOOT, CAT.7 SIZE 24CM SL	F20-S3-07A24-SL
722402472	SHOCKWAVE FOOT, CAT.7 SIZE 24CM SR	F20-S3-07A24-SR
722402570	SHOCKWAVE FOOT, CAT.7 SIZE 25CM RU	F20-S3-07A25-RU
722402571	SHOCKWAVE FOOT, CAT.7 SIZE 25CM SL	F20-S3-07A25-SL
722402572	SHOCKWAVE FOOT, CAT.7 SIZE 25CM SR	F20-S3-07A25-SR
722402670	SHOCKWAVE FOOT, CAT.7 SIZE 26CM RU	F20-S3-07A26-RU
722402671	SHOCKWAVE FOOT, CAT.7 SIZE 26CM SL	F20-S3-07A26-SL
722402672	SHOCKWAVE FOOT, CAT.7 SIZE 26CM SR	F20-S3-07A26-SR
722402770	SHOCKWAVE FOOT, CAT.7 SIZE 27CM RU	F20-S3-07A27-RU
722402771	SHOCKWAVE FOOT, CAT.7 SIZE 27CM SL	F20-S3-07A27-SL
722402772	SHOCKWAVE FOOT, CAT.7 SIZE 27CM SR	F20-S3-07A27-SR

722402870	SHOCKWAVE FOOT, CAT.7 SIZE 28CM RU	F20-S3-07A28-RU
722402871	SHOCKWAVE FOOT, CAT.7 SIZE 28CM SL	F20-S3-07A28-SL
722402872	SHOCKWAVE FOOT, CAT.7 SIZE 28CM SR	F20-S3-07A28-SR
722402970	SHOCKWAVE FOOT, CAT.7 SIZE 29CM RU	F20-S3-07A29-RU
722403070	SHOCKWAVE FOOT, CAT.7 SIZE 30CM RU	F20-S3-07A30-RU
722402380	SHOCKWAVE FOOT, CAT.8 SIZE 23CM RU	F20-S3-08A23-RU
722402381	SHOCKWAVE FOOT, CAT.8 SIZE 23CM SL	F20-S3-08A23-SL
722402382	SHOCKWAVE FOOT, CAT.8 SIZE 23CM SR	F20-S3-08A23-SR
722402480	SHOCKWAVE FOOT, CAT.8 SIZE 24CM RU	F20-S3-08A24-RU
722402481	SHOCKWAVE FOOT, CAT.8 SIZE 24CM SL	F20-S3-08A24-SL
722402482	SHOCKWAVE FOOT, CAT.8 SIZE 24CM SR	F20-S3-08A24-SR
722402580	SHOCKWAVE FOOT, CAT.8 SIZE 25CM RU	F20-S3-08A25-RU
722402581	SHOCKWAVE FOOT, CAT.8 SIZE 25CM SL	F20-S3-08A25-SL
722402582	SHOCKWAVE FOOT, CAT.8 SIZE 25CM SR	F20-S3-08A25-SR
722402680	SHOCKWAVE FOOT, CAT.8 SIZE 26CM RU	F20-S3-08A26-RU
722402681	SHOCKWAVE FOOT, CAT.8 SIZE 26CM SL	F20-S3-08A26-SL
722402682	SHOCKWAVE FOOT, CAT.8 SIZE 26CM SR	F20-S3-08A26-SR
722402780	SHOCKWAVE FOOT, CAT.8 SIZE 27CM RU	F20-S3-08A27-RU
722402781	SHOCKWAVE FOOT, CAT.8 SIZE 27CM SL	F20-S3-08A27-SL
722402782	SHOCKWAVE FOOT, CAT.8 SIZE 27CM SR	F20-S3-08A27-SR
722402880	SHOCKWAVE FOOT, CAT.8 SIZE 28CM RU	F20-S3-08A28-RU
722402881	SHOCKWAVE FOOT, CAT.8 SIZE 28CM SL	F20-S3-08A28-SL
722402882	SHOCKWAVE FOOT, CAT.8 SIZE 28CM SR	F20-S3-08A28-SR
722402980	SHOCKWAVE FOOT, CAT.8 SIZE 29CM RU	F20-S3-08A29-RU
722403080	SHOCKWAVE FOOT, CAT.8 SIZE 30CM RU	F20-S3-08A30-RU
