


HOPPER



The HOPPER running blade is designed for running on all types of terrain and at all distances by people of all levels, athletes and beginners alike. Its comfort, versatility and ease of use mean it's quick to get to grips with a wide range of activities.

Two blade sizes available: AUDACIEUSE (>24cm Build height) and INTREPID (>19cm build height).

- Versatility and dynamism
- Comfort and flexibility
- Eco-responsibility

TECHNICAL FEATURES

Model	Audacieuse	Intrepide
Item N°.	HL01	HL02
Weight	840 g	740 g
Max. patient weight	45 / 110 kg	35 / 110 kg
Minimal clearance	240 mm	190 mm
Material	Recycled carbon fiber from the aerospace industry	
Outsole	Salomon® Blue Fire outsole	



HS01B



HA01F

HA01M


1 = Low 2 = Medium 3 = High

HL	Category
Model	See table
Audacieuse : 01	
Intrepid : 02	

Category selection according to weight and activity													
Weight		77-100 lbs 35-45 kg		101-121 lbs 45-55 kg		122-148 lbs 55-67 kg		149-172 lbs 68-78 kg		173-199 lbs 79-90 kg		200-242 lbs 90-110 kg	
Activity	Minimal Clearance*	1/2	3	1/2	3	1/2	3	1/2	3	1/2	3	1/2	3
Audacieuse	240 mm			XS	S	S	M	M	L	L	XL	XL	XL
Intrepide	190 mm	PXS	PXS	PXS	PS	PS	PM	PM	PL	PL	PXL	PXL	PXL

For bilateral patients: increase by one category for added stability.

*Minimal clearance: distance between the bottom of the socket and the ground on the walking prosthesis with shoe.



1X998 To avoid damage to the pyramid during trials

